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
Connecting LAN cables

Connect the LAN cable to the LAN connector on the underside of the device.

Connect the other end of the LAN cable to the hub to which the receiver PC is connected.

CAUTION

- Be sure to connect to the PC via a network hub.
Do not connect to the PC directly. Connection failure may result.

2.4 Confirmation and Setting of Device

2.4.1 Setting the Device Transmission

Set the network parameters of the device based on the information of “2.3.3 Confirming Items Necessary for Other Settings” (page 44).

1 Turn on the device to enter the PARAMETER SETTING screen.

2 Select “NETWORK1” from the parameter items.
The NETWORK1 setting screen appears.

3 Set the “111. NETWORK” parameter to “YES”.
For connection to NAVIS, set “ACK”.

PARAMETER	SETTING
[NETWORK1]	
111. NETWORK	YES
112. DHCP	NO
113. IP	192. 168. 0. 30
114. MASK	255. 255. 255. 0
115. USER	GUEST
116. PASSWORD	
117. DOMAIN	WORKGROUP
118. PC NAME	PC
119. FOLDER	DATA

4 Set the “112. DHCP” parameter to “YES” or “NO”.

For cases where an IP address is automatically assigned, use the DHCP parameter (select “YES”).

When the “112. DHCP” parameter is set to “YES”, Steps 5 to 10 are not necessary to be conducted. Proceed to Step 11.

When an IP address is issued, select “NO”.

5 When the “113. IP” parameter is selected, the enter (✓) button is displayed.

PARAMETER	SETTING
[NETWORK1]	
111. NETWORK	YES
112. DHCP	NO
113. IP	192. 168. 0. 30
114. MASK	255. 255. 255. 0
115. USER	GUEST
116. PASSWORD	
117. DOMAIN	WORKGROUP
118. PC NAME	PC
119. FOLDER	DATA

6 Press the enter (✓) button to display the IP address setting screen.

IP ADDRESS:
192. 168. 0. 30


CHANGE SET

- 7** Input the IP address based on the information of “2.3.3 Confirming Items Necessary for Other Settings” (page 44) by the following button operation.

(ex.: 10. 50. 1. 10)

 Forward movement /  Back-ward movement or pressing the start button	Moves the cursor to the right or left. Pressing the start button moves the cursor to the right.
Rotating the upper part of the joystick	Increases/decreases the numeric value.

- 8** Pressing the exit  button displays a screen to confirm the input IP address.

- 9** Select “YES” and press the enter  button.


Select “YES” or “NO” by pressing the  /  button or rotating the joystick.

The screen returns to the “NETWORK 1” setting screen.

- 10** Select the “114. MASK” parameter and set the subnet mask in the same manner as Steps 5 to 9.


(ex.: 255. 255. 255. 0)

- 11** Select the “115. USER” parameter.

- 12** When the “115. IP” parameter is selected, the enter  button is displayed.

- 13** Press the enter  button to display the user setting screen.

- 14** Input the login user name based on the information of “2.3.3 Confirming Items Necessary for Other Settings” (page 44) by the following button operation.

* Entry is case insensitive. Input all with uppercase letters.

(ex.: NIDEK)

 Forward movement /  Backward movement	Moves the character to be input (cursor) in the character list field to the right or left.
Rotating the upper part of the joystick	Moves the cursor in the input field to the right or left.
Pressing the start button	Determines the character to be input and moves the cursor to the next position in the character list field.
 Clear	Erases the character at the cursor position in the character list field.

- 15** Pressing the exit  button displays a screen to confirm the input user name.

- 16** Select “YES” and press the enter  button.

The screen returns to the “NETWORK 1” setting screen.

- 17** Select the “116. PASSWORD” parameter and then set the login password name in the same manner as Steps 11 to 16 based on the information of “2.3.3 Confirming Items Necessary for Other Settings” (page 44).

* Entry is case sensitive. Input uppercase and lowercase letters correctly.

(ex.: NIDEK)

- 18** Select the “117. DOMAIN” parameter and then set the domain name in the same manner as Steps 11 to 16 based on the information of “2.3.3 Confirming Items Necessary for Other Settings” (page 44).

(ex.: NAVIS)

(If no domain exists, confirm the work group name.) (ex.: WORKGROUP)


- 19** Select the “118. PC NAME” parameter and then set the PC name in the same manner as Steps 11 to 16 based on the information of “2.3.3 Confirming Items Necessary for Other Settings” (page 44).

(ex.: 10. 50. 1. 1)

- 20** Select the “119. DOMAIN” parameter and then set the shared folder of the PC to which measurement data is sent in the same manner as Steps 11 to 16 based on the information of “2.3.3 Confirming Items Necessary for Other Settings” (page 44).

(ex.: ARK)

- 21** Select “NETWORK2” from the parameter items and then set the “121. RING IMAGE”, “122. RETRO IMAGE”, and “123. ACC GRAPH” parameters for transmission of image data.

- 22** After all “NETWORK1” to “NETWORK2” parameters are set, press the exit  button twice to exit from the PARAMETERSETTING screen.

- Confirm that the PARAMETER SETTING screen is completely exited and the normal measurement screen is displayed.
- Unless the PARAMETER SETTING screen is exited, network settings are not saved.

- 23** Turn off the device and then turn it on again to reboot.

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