Geräteanbindung OCULUS/NIDEK LM-1000P

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1	Basic specifications	In accordance with the RS-232C specifications
2	Communication method	Asynchronous
3	Transmission type	Half duplex
4	Communication mode (ComMode)	<u>NIDEK</u> , PC NCP10
5	Baud rate (BaudRate)	1200, 2400, 4800, <u>9600</u> (bps)
6	Parity (Parity check)	OFF, <u>ODD</u> , EVEN
7	Data bits	7bit, <u>8bit</u>
8	Stop bits	<u>1bit</u> , 2bit
9	CR code	<u>OFF</u> ON
10	Data code	ASCII code

The underlined options show the factory settings for communication with NIDEK devices.

The settings of Parameter 4 ComMode to Parameter 9 CR Code are changeable.

For the procedure for changing the parameter settings, see "2.13 Setting Parameters" of the LM-1000, LM-1000P, or LM-1200 Operator's Manual.

The settings of Parameter 5 BaudRate to Parameter 8 StopBits are changeable only when Parameter 4 ComMode is set to PC or NCP10.

If Parameter 4 ComMode is set to NIDEK, Parameter 5 BaudRate to Parameter 8 StopBits are automatically set for connection with NIDEK devices; once the parameters are set automatically, the settings are unchangeable.

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- Addition of CR code

Set the CR Code parameter to ON when the communication software on the computer side needs the CR code to terminate the reception of the character strings. (BASIC language etc.)

The CR code is added to the data sent from the LM; For character strings sent from the computer, whether the CR code is added or not does not matter.

- To enable communication with a computer, it is necessary to keep consistency in the settings of communication parameters between the LM and computer.

Before data communication, also confirm the settings of the communication parameters on the computer side.

○ Setting related to the transmission of prism data

Set the Prism Tx parameter to suit the needs for prism data.


11	Prism Tx	OFF, ON, <u>Display</u>
----	----------	-------------------------

OFF	The measured prism data is not transmitted at any time. Choose this option not to transmit prism data or restrain the size of data to be transmitted.
ON	The measured prism data is transmitted at all times. Choose this option to make the external device control and align measured lenses.
Display	The measured prism data is transmitted in a coordinate system in which the data is represented on the screen, OFF: The measured prism data is not transmitted. BU/D BI/O: The measured prism data in a rectangular coordinates is transmitted. P-B: The measured prism data in a polar coordinates is transmitted.

- Set the representation form of prism data by setting the Prism parameter.

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Connect signal wires between the LM and the external device as shown above.

Be sure to connect a wire between RTS (7) and CTS (8) on the computer side.

It is unnecessary connect a wire between RTS (4) and CTS (5) on the LM side.

Use the DIN 8-pin connector (DIN Standard No.45329) on the LM side.

In the above example, the D-Sub 9-pin connector is connected to the computer side.


If the used connector is different from the D-Sub 9-pin connector, connect respective signal wires as shown by referring to the instruction manual provided with the computer.

Use a special interface cable that is available as an option (OPIF-11, 5 m) between the DIN 8-pin and D-Sub 9 pin connectors if necessary.

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
The screenshot shows the 'Patienten-ID / Kennungen' tab. Fields include:

- Exportformat der CGM MEDISTAR Patientennummer (GDT-Kennung 3000)
- Importformat der CGM MEDISTAR Patientennummer (GDT-Kennung 3000)
- Importformat führende Nullen unterdrücken (checkbox)
- Gesendete GDT-Version (GDT-Kennung 9218): 02.10
- Gesendete GDT-Sender-Kennung (GDT-Kennung 8316): EDV1
- Gesendete GDT-Empfänger-Kennung (GDT-Kennung 8315): [Redacted]

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Untersuchungsart OPT003, alternativ **LENSM01** (nicht GDT-Standard).

The screenshot shows the 'Export' tab. Settings include:

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The screenshot shows the 'MD-Eintrag' tab. Settings include:

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